

Data Sheet

CellGenix® Recombinant Human Interleukin-1 beta (rh IL-1β)

Preclinical Grade – Order No.: 1411-010 (10 µg)*, 1411-050 (50 µg)

Product Characteristics

Source	<i>E. coli</i>
Description	Human Interleukin-1 beta, accession # P01584, Ala117-Ser269 N-terminal Met Molecular mass 17.5 kDa
Formulation	Lyophilized from a 0.2 µm-filtered solution containing 1.5 mM potassium phosphate, 8.1 mM sodium phosphate, 2.7 mM potassium chloride, and 137 mM sodium chloride, pH 7.5.
Intended use	For preclinical <i>ex vivo</i> use. Not intended for therapeutic use.

Quality Parameters

Activity	≥ 100 x 10 ⁶ IU/mg calibrated against NIBSC #86/680 Measured in a cell proliferation assay using an IL-1β-dependent cell line, RPMI-1788
Purity	≥ 95 %, as determined by SDS-PAGE (under reducing and non-reducing conditions, visualized by silver staining)
Endotoxin	< 1000 EU/mg, as determined by LAL gel clot test
Sterility	Sterility test of the vial product
Mass per vial	1411-010: 10 µg, 1411-050: 50 µg
Animal-derived component-free	ADCF Level 2: The final product contains neither animal- nor human-derived materials. ADCF Level 2 cytokines are produced in our dedicated animal-free facility. No animal-derived components are used throughout the complete production process. All ADCF Level 2 cytokines are produced in <i>E. coli</i> .

Shipment & Storage

Transport	Ambient temperature. Please refer to Technote (www.cellgenix.com)
Shelf life	3 years from date of shipment
Storage & Stability	Store lyophilized cytokine at -20 °C to -80 °C. Avoid repeated freeze/thaw cycles.

Handling Instructions

Reconstitution	Recommended in sterile water to a final concentration of 100 µg/ml (for 10 µg vials) or 250 µg/ml (for 50 µg vials).
Dilution	Recommended in CellGenix® serum-free media. For dilution with protein free medium, a carrier protein (0.1-1 % albumin or 1-10 % appropriate serum) has to be included. Failure to dilute product according to these instructions may result in loss of activity.

Quality Statement

Final manufacturing steps and QC are performed in a GMP facility. No animal- or human-derived components are present in the final product and no animal- or human-derived materials were used in production (ADCF Level 2).